




### Specifications

Electrical requirements	110VAC 60Hz single phase, 1200 watts
Chamber dimensions	14 x 12 x 14" (WxDxH)
Unit dimensions	20 x 23 x 26.5" (WxDxH)
Shipping dimensions	27 x 24 x 31" (WxDxH)
Temperature range	Ambient temperature to 300°C
Temperature uniformity	+/- 2°F
Temperature controller precision	0.1°F
Heating elements	In the back of the oven
Timer range	1 to 9999 minutes
Racks included	2 (additional racks available at extra cost)
Maximum rack capacity	5 (height between each rack is 2 inches)
Exhaust port	2" diameter x 1, adjustable
Shipping weight	110 Lb
Access port	2" diameter with cover, one
Certification	CE, UL/CSA (optional)

### Standard package

Part description	Quantity	Part image
AI 1.5 cu ft forced air oven	1 pc	
Stainless steel rack and brackets	2 sets	
110 volt 3-prong power cord	1 pc	
User's manual	1 pc	