

## RAPGEF5 Polyclonal Antibody

Catalog number: 12556-1-AP


Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RAPGEF5; GFR; KIAA0277; MR-  
GEF; REPAC



human brain tissue were subjected to SDS PAGE followed by western blot with 12556-1-AP(RAPGEF5 Antibody) at dilution of 1:600

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse; other species not tested.
Positive WB detected in:	Human brain tissue, mouse brain tissue
Calculated RAPGEF5 MW:	580aa,67.7kd
Observed RAPGEF5 MW:	68-75 kDa
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3229
GenBank accession number:	BC039203
Gene ID (NCBI):	9771
Full name:	Rap guanine nucleotide exchange factor (GEF) 5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.