

Product Name:	ABS SSC 47-6
Description:	<p> ABS SSC Class Includes 6 Canisters Total 2ml Vial Capacity: 1320 Liquid Nitrogen Capacity: 47.4 Liters ½ CC Straw (1 Level Bulk): 6200 ½ CC Straw (10 Cane): 3900+ Static Evaporation Rate (Liters/day): 0.39 Normal Working Time: 76 days Neck Diameter (mm): 127 Exterior Height (mm): 673 Exterior Diameter (mm): 508 Canister Diameter (mm): 101.6 Canister Height (mm): 279.4 Durable aluminum tank with powder coated exterior Lid lock provision fiberglass restricted opening neck design Compatible design for both liquid and vapor phase storage LN2 measuring rod for level validation Color-coded canisters for inventory identification One-year parts and labor warranty Three-year vacuum warranty </p>