

CLK2 Polyclonal Antibody

Catalog number: 11092-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CLK2; MGC61500; hCLK2

Applications

Tested applications: ELISA
Species specificity: Human, mouse, rat; other species not tested.
Calculated CLK2 MW: 60kd
Observed CLK2 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag1549
GenBank accession number: BC014067
Gene ID (NCBI): 1196
Full name: CDC-like kinase 2

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.