

KCNIP4 Polyclonal Antibody

Catalog number: 13748-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

KCNIP4; CALP; KCHIP4;

MGC44947



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13748-1-AP(KCNIP4 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse brain tissue
Calculated KCNIP4 MW:	250aa,28kd
Observed KCNIP4 MW:	29kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4688
GenBank accession number:	BC032520
Gene ID (NCBI):	80333
Full name:	Kv channel interacting protein 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.