

LabStrong Replacement Cartridges & Parts for Thermo Scientific™ Barnstead™ Water Systems


REDESIGNED!

High Availability

Low Prices




UV Lamps



Cartridges



RO Membranes



O-Rings



Filters

Exceptional Quality



Heating Elements



Pumps

Maintaining lab equipment is vital to lab productivity. Keep equipment running smoothly with quality products!

Wide Variety

- Suitable for B-Pure™, Bantam™, E-Pure™, EASYpure™, Diamond™ RO, Infinity, Mega-Pure™, NANOpure™ Diamond, NANOpure II, NANOpure Quad, TII, ROpure ST and many other water systems
- Hose nipple cartridges available

Budget Friendly

- Outstanding quality offered at economical price

High Quality

- Comparable to the Thermo Scientific Barnstead line of products
- Many items feature new and improved design

Simple to Order

- Sold on individual basis to fill any need
- Most items in stock and can ship same day at no additional charge

 - MADE IN USA OF USA AND IMPORTED PARTS - 

Contact Us

Phone: (563) 588-8900 • Fax: (563) 588-3601
sales@labstrong.com • www.labstrong.com

 **LabStrong™**
Strength Through Innovation

LabStrong #	Barnstead #	Description	Where Used	List Price
D0803LS	D0803	High Capacity (Two-Bed)	B-Pure, Bantam, Cooling Loops, E-Pure, NANOpure, NANOpure II & NANOpure Quad	\$131.00
D0809LS	D0809	Ultrapure (Mixed-Bed)	B-Pure, Bantam, Cooling Loops, E-Pure, NANOpure, NANOpure II & NANOpure Quad	\$134.00
D0813LS	D0813	Organic Removal (Carbon)	B-Pure, Bantam & Cooling Loops	\$115.00
D0815LS	D0815	Cation Removal (Cation) NEW!	B-Pure, Bantam & Cooling Loops	\$109.00
D0820LS	D0820	Organic Free	NANOpure II & older NANOpure	\$164.00
D0832LS	D0832	Pretreatment Cartridge (Mixed-Bed/Carbon)	B-Pure & Still Pretreatment	\$120.00
D0835LS	D0835	Combination Cartridge (Colloid/Organic)	E-Pure, older NANOpure, NANOpure II & NANOpure Quad	\$134.00
D0836LS	D0836	Macropure (Macroreticular/Carbon)	E-Pure, older NANOpure, NANOpure II & NANOpure Quad	\$175.00
D8901LS	D8901	HN (3/8") Dyed High Capacity (Two-Bed)	Hose Nipple, non-pressure application	\$134.00
D8902LS	D8902	HN (3/8") Dyed Ultrapure (Mixed-Bed)	Hose Nipple, non-pressure application	\$134.00
D8904LS	D8904	HN (3/8") Organic Removal - No Dye (Carbon)	Hose Nipple, non-pressure application	\$107.00
D8905LS	D8905	HN (3/8") Dyed Cation Removal NEW!	Hose Nipple, non-pressure application	\$107.00
D8911LS	D8911	HN (3/8") Ultrapure - No Dye (Mixed-Bed)	Hose Nipple, non-pressure application	\$130.00
D8921LS	D8921	HN (3/8") Dyed Pretreatment (Mixed-Bed/Carbon) NEW!	Hose Nipple, non-pressure application	\$130.00
D8922LS	D8922	HN (3/8") Dyed Organic & Ultrapure Removal (Mixed-Bed/Carbon)	Hose Nipple, non-pressure application	\$140.00
D8950LS	D8950	HN (3/8") 2/3 Size Dyed High Capacity (Two-Bed) NEW!	Hose Nipple, non-pressure application	\$113.00
D8951LS	D8951	HN (3/8") 2/3 Size Dyed Organic Removal (Mixed-Bed/Carbon) NEW!	Hose Nipple, non-pressure application	\$120.00
D5021LS	D5021	Organic Free	E-Pure & NANOpure Quad	\$171.00
D5027LS	D5027	Ultrapure (Mixed-Bed)	E-Pure, NANOpure Quad & ROpure ST	\$133.00
D5028LS	D5028	E-Pure Cartridge Kit (Tap Feed/4 Mod Standard)	E-Pure	\$560.00
D5023LS	D5023	E-Pure Cartridge Kit (Tap Feed/4 Mod Organic Free)	E-Pure	\$642.00
D50227LS	D50227	E-Pure Cartridge Kit (Pretreated Feed/4 Mod Standard)	E-Pure	\$561.00
D50228LS	D50228	E-Pure Cartridge Kit (Pretreated Feed/4 Mod Organic Free)	E-Pure	\$642.00
D5029LS	D5029	E-Pure Cartridge Kit (3 Mod Standard)	E-Pure	\$434.00
D5022LS	D5022	E-Pure Cartridge Kit (3 Mod Organic Free)	E-Pure	\$500.00
D4801LS	D4801	NANOpure Quad Style Cartridge Kit (Tap Feed/Standard)	NANOpure Quad	\$561.00
D4802LS	D4802	NANOpure Quad Style Cartridge Kit (Tap Feed/Organic Free)	NANOpure Quad	\$658.00
D5025LS	D5025	NANOpure Quad Style Cartridge Kit (Pretreated Feed/Standard)	NANOpure Quad	\$561.00
D5026LS	D5026	NANOpure Quad Style Cartridge Kit (Pretreated Feed/Organic Free)	NANOpure Quad	\$704.00
D3804LS	D3804	NANOpure II Cartridge Kit (4 Mod Standard)	NANOpure II & Older NANOpure	\$589.00
D3805LS	D3805	NANOpure II Cartridge Kit (3 Mod Standard)	NANOpure II & Older NANOpure	\$448.00
D3806LS	D3806	NANOpure II Cartridge Kit (3 Mod Organic Free)	NANOpure II & Older NANOpure	\$499.00
D3807LS	D3807	NANOpure II Cartridge Kit (4 Mod Organic Free)	NANOpure II & Older NANOpure	\$652.00
D400377LS	D400377	Mega-Pure High Purity Cartridge NEW!	Mega-Pure D1, D2 & MP12-A Stills	\$256.00
D400499LS	D400499	Mega-Pure High Capacity Cartridge NEW!	Mega-Pure D1, D2 & MP12-A Stills	\$268.00
D440265LS	D440265	Mega-Pure Organic and Chlorine Removal Cartridge NEW!	Mega-Pure D2 & MP12-A Stills	\$235.00
LabStrong #	Barnstead #	Description	Where Used	List Price
06808LS	06808 & 06808-DB	NANOpure II Canister O-Ring NEW!	NANOpure II	\$6.00
25001LS	25001 & 25001-DB	Air Filter	Water Storage Reservoirs (Cabinetized Stills, Classic Stills, Diamond TII, etc.)	\$65.00
410671LS	410671	Mega-Pure Gasket 4 Pack NEW!	Mega-Pure Single (D1) & Dual (D2) Deionizer Systems & Mega-Pure MP-12A Stills	\$70.00
740880LS	740880	Mega-Pure MP-1 Heating Element NEW!	Mega-Pure Stills (MP-1)	\$620.00
740883LS	740883 & 740895	Mega-Pure Heating Element	Mega-Pure Stills (MP-3A, MP-6A, MP-11A & MP-12A)	\$620.00
CMX25LS	CMX25	Cleaning Syringe NEW!	EASYpure II (LF, UV & UV/UF), NANOpure Infinity & NANOpure Diamond	\$35.00
D8900LS	D8900	Wall Mount Bracket	Hose Nipple Cartridges	\$42.00
FL0016-1	D3750	Sterile 0.2 Micron Final Filter with Endbell	B-Pure, E-Pure, all NANOpure systems (including Diamond), AccuDispense, EASYpure, all EASYpure II systems, EASYpure RoDi & Diamond TII	\$99.00
FL1265X1LS	FL1265X1	Encapsulated RO Membrane	Diamond RO & Diamond TII	\$472.00
GSX27LS	GSX27	Cartridge O-Ring	B-Pure, E-Pure & NANOpure Quad	\$4.00
GSX28LS	GSX28	Canister O-Ring	B-Pure, E-Pure & NANOpure Quad	\$6.00
LMP13	LMX13	185/254 nm UV Lamp	NANOpure Quad, Infinity, EASYpure, EASYpure II (RF/UV, UV & UV/UF), EASYpure RoDi, older NANOpure Diamond & Diamond TII	\$142.00
LMP01	LMX13	185/254 nm UV Lamp	Diamond TII & NANOpure Diamond with new style ballast	\$142.00
PU1190X1LS	PU1190X1	Diaphragm Pump	NANOpure Infinity Retrofit Kit (PU896X3) & NANOpure Diamond	\$427.00