

- For Emergency Use Authorization (EUA) only.
- This test can be used at home on people aged 2 years old and up.
- Items necessary to use the kit, but not provided:
 - Timer
- For Symbol Glossary, refer to Instructions for Use.

*In the USA, this product has not been FDA cleared or approved; but has been authorized by FDA under an Emergency Use Authorization. This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of *in vitro* diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.



Insert tube here
Inserte el tubo aquí



Insert tube here
Inserte el tubo aquí



If you have symptoms of COVID-19, you can use a single test.
If you do not have symptoms of COVID-19, you will need at least two tests per person.
You may need to purchase additional tests to perform serial (repeat) testing.
This test is more likely to give you a false negative result when you have COVID-19 than a lab-based molecular test.

The CLINITEST® Rapid COVID-19 Antigen Self-Test is an immunochromatographic membrane assay that uses antibodies to detect nucleocapsid protein from SARS-CoV-2 in anterior nasal swabs from those with symptoms of COVID-19 within the first 7 days of symptom onset or those with or without symptoms or other epidemiological reasons to suspect COVID-19 when tested twice over three days with at least 24 hours (and no more than 48 hours) between tests.

LOT



REF GCCOV-502a-H5US
(11556712)



00816490025740



***.clinitest.siemens-healthineers.com/us

B01601-01

SIEMENS
Healthineers

CLINITEST®

Rapid COVID-19 Antigen Self-Test

Contents of the kit:
5 Test Devices
5 Sterile Swabs
5 Extraction Tubes with Buffer and Tips
Instructions for Use



Distributed by Siemens Healthineers



Healgen Scientific Limited Liability Company
Address: 3818 Fuqua Street, Houston, TX 77047, USA.
Tel: +1 713-733-8088 Fax: +1 713-733-8848
Website: www.healgen.com