

SMSLC SERIES


SMS (SPUNBOND MELTBLOWN SPUNBOND) LAB COAT

Introducing our SMS Lab Coat - Our SMS Lab Coat offers top-tier protection in labs, clean rooms, and similar environments. Designed with your safety in mind, our SMS Lab Coat offers excellent protection. Made from high-quality SMS (Spunbond Meltblown Spunbond) fabric, it effectively shields you from particles and fluids. The lightweight yet durable material ensures comfort, while the button-up front and elastic wrists provide a secure fit. Our SMS Lab Coat is perfect for both men and women, making it a versatile choice for various work environments. Prioritize safety without sacrificing comfort. Choose our SMS Lab Coat for reliable protection in the lab and clean room.



SIZE	PACKAGING	THOMAS NO. SKU
Small	10pcs/polybag,30pcs/case	24A00G494 SMSLC101
Medium	10pcs/polybag,30pcs/case	24A00G495 SMSLC102
Large	10pcs/polybag,30pcs/case	24A00G496 SMSLC103
X-Large	10pcs/polybag,30pcs/case	24A00G497 SMSLC104
2XL	10pcs/polybag,30pcs/case	24A00G498 SMSLC105
3XL	10pcs/polybag,30pcs/case	24A00G499 SMSLC106

in inches



Size	A. Length	B. Chest	C. Shoulder	D. Sleeve	E. Cuff		F. Distance Between Snaps	H. Sleeve Width
					Relax	Stretch		
S	37.1	24.0	7.6	22.2	3.4	5.5	5.9	10.0
M	38.3	24.8	7.8	22.5	3.5	6.0	5.9	10.4
L	39.4	25.5	8.1	23.0	3.5	6.0	5.9	10.8
XL	40.6	26.5	8.3	23.5	3.5	6.0	5.9	10.8
2XL	41.8	28.0	8.7	24.5	3.5	6.0	5.9	11.2
3XL	43.0	29.5	8.9	25.2	3.8	6.5	5.9	11.6
4XL	44.2	30.0	9.1	26.0	4.0	7.5	5.9	12.0
5XL	45.4	31.0	9.3	26.5	4.3	8.0	5.9	12.4
Deviation:	1.00+/-	1.00+/-	0.50+/-	1.00+/-	0.50+/-	0.50+/-	0.50+/-	0.50+/-