

TBRG4 Polyclonal Antibody

Catalog number: 10177-2-AP

Size: 93 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TBRG4; CPR2, FASTKD4,

KIAA0948, Protein TBRG4,

TBRG4



A431 cells were subjected to SDS PAGE followed by western blot with 10177-2-AP (TBRG4 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated TBRG4 MW:	71 kDa
Observed TBRG4 MW:	58 kDa
Positive WB detected in	A431 cells, HepG2 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0212
GenBank accession number:	BC002732
Gene ID (NCBI):	9238
Full name:	Transforming growth factor beta regulator 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.