

For Research Use Only

INTS12 POLYCLONAL ANTIBODY



Catalog Number: 16455-1-AP

Basic Information

Catalog Number:
16455-1-AP

Size:
21 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification

Immunogen Catalog Number:
AG9403

GenBank Accession Number:
BC014442

GeneID (NCBI):
57117

Full Name:
integrator complex subunit 12

Calculated MW:
462aa, 49 kDa

Observed MW:
49-52kd

Recommended Dilutions:
WB 1:1000-1:10000

Applications

Tested Applications:
ELISA, WB

Species Specificity:
human, mouse, rat

Positive Controls:

WB : HepG2 cells;

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
--------	-----------	---------	-------------

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 16455-1-AP(INTS12 antibody) at dilution of 1:1000