

CaMKIV Polyclonal Antibody

Catalog number: 13263-1-AP

Size: 45 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

CAMK4; CaM kinase GR, CaMK

GR, CAMK4, CaMKIV



human cerebellum tissue were subjected to SDS PAGE followed by western blot with 13263-1-AP(CAMK4 antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated CaMKIV MW:	473aa,52kd
Observed CaMKIV MW:	65kd
Positive WB detected in	Human cerebellum tissue, mouse thymus tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3922
GenBank accession number:	BC025687
Gene ID (NCBI):	814
Full name:	Calcium/calmodulin-dependent protein kinase IV

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.