

TPRKB Polyclonal Antibody

Catalog number: 12800-1-AP


Size: 40 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TPRKB; CGI-121



human heart tissue were subjected to SDS PAGE followed by western blot with 12800-1-AP (TPRKB Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human heart tissue
Calculated TPRKB MW:	175aa, 20kd
Observed TPRKB MW:	20kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3551
GenBank accession number:	BC029492
Gene ID (NCBI):	51002
Full name:	TP53RK binding protein

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.