

ATP1A3 Polyclonal Antibody

Catalog number: 10868-1-AP

Size: 50 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ATP1A3; ATP1A3, DYT12, RDP,

Sodium pump subunit alpha 3

Applications

Tested applications: ELISA
Species specificity: Human, Mouse, Rat; other species not tested.
Calculated ATP1A3 MW: 113 kDa
Observed ATP1A3 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Ag1313
GenBank accession number: BC015566
Gene ID (NCBI): 478
Full name: ATPase, Na⁺/K⁺ transporting, alpha 3 polypeptide

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.