

FNTA Polyclonal Antibody

Catalog number: 12274-1-AP

Size: 29 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

FNTA; FNTA, FPTA, FTase alpha,

GGTase I alpha, PGGT1A,

PTAR2



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 12274-1-AP (FNTA antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated FNTA MW:	44 kDa
Observed FNTA MW:	44 kDa
Positive WB detected in	NIH/3T3 cells, 3T3-L1 cells, A431 cells, Jurkat cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2918
GenBank accession number:	BC017029
Gene ID (NCBI):	2339
Full name:	Farnesyltransferase, CAAX box, alpha

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.