

TBL1X Polyclonal Antibody

Catalog number: 13540-1-AP


Size: 24 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TBL1X; EBI; SMAP55; TBL1



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with 13540-1-AP(TBL1X Antibody) at dilution of 1:1200

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse pancreas tissue
Calculated TBL1X MW:	526aa,57kd
Observed TBL1X MW:	57kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4477
GenBank accession number:	BC032708
Gene ID (NCBI):	6907
Full name:	Transducin (beta)-like 1X-linked

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.