

GAN Polyclonal Antibody

Catalog number: 14305-1-AP


Size: 42 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GAN; FLJ38059; GAN1; KLHL16



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14305-1-AP(GAN Antibody) at dilution of 1:1500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse brain tissue, HEK-293 cells
Calculated GAN MW:	68kd
Observed GAN MW:	68kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag5421
GenBank accession number:	BC044840
Gene ID (NCBI):	8139
Full name:	Gigaxonin

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.