

B4GALT7 Polyclonal Antibody

Catalog number: 10535-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

B4GALT7; B4GAL-T7; XGALT-1;

XGALT1; XGPT1; beta4Gal-T7



MCF7 cells were subjected to SDS PAGE followed by western blot with 10535-1-AP (B4GALT7 Antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Cited applications:	WB
Species specificity:	Human; other species not tested.
Cited species:	Human
Positive WB detected in:	MCF7 cells, HeLa cells
Calculated B4GALT7 MW:	37kd
Observed B4GALT7 MW:	37kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0793
GenBank accession number:	BC007317
Gene ID (NCBI):	11285
Full name:	Xylosylprotein beta 1,4-galactosyltransferase, polypeptide 7 (galactosyltransferase I)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.