

CDK2AP1 Polyclonal Antibody

Catalog number: 13060-2-AP

Size: 25 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

CDK2AP1; CDK2 associated

protein 1, CDK2AP1, CDKAP1,

Deleted in oral cancer 1, doc 1,

DOC1, DORC1, p12DOC 1, ST19



HEK-293 cells were subjected to SDS PAGE followed by western blot with 13060-2-AP (CDK2AP1 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated CDK2AP1 MW:	115aa, 13 kDa
Observed CDK2AP1 MW:	12 kDa
Positive WB detected in	HEK-293 cells, HeLa cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3736
GenBank accession number:	BC034717
Gene ID (NCBI):	8099
Full name:	Cyclin-dependent kinase 2 associated protein 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.