







2-Stage Horizontal Rotor

Our patented swing-out horizontal rotor design incorporates a unique test tube holder that produces the ideal horizontally separated samples.

Loading Angle



Rotation Angle



DRUCKER MODEL 642E SINGLE-SPEED ELECTRONIC CENTRIFUGE

The 642E is the **world's best-selling clinical laboratory centrifuge**. Easy two-button operation ensures repeatability and standardization, while advanced design delivers high-quality blood separation every time. Its preset speed and runtime make the 642E the ideal choice for consistency in doctor's offices, patient service centers, and outreach locations.

Easy to Use

Spin samples with the press of just one button. The 642E processes 75-100 mm (3-10 mL) tubes for chemistry, coag, and PPP.

Simple Settings

The 642E features a preset runtime of 10 minutes at 1,600 xg for all common chemistry needs.

Built to Last

No routine maintenance. The brushless motor and unique convection-based cooling system ensure years of high-quality separations.

Engineered for Safety

Protect your lab with the clear, shatterproof lid, which locks until the centrifuge has stopped completely.

A Great Value

Get the world's #1 clinical laboratory centrifuge at a price that fits any budget. The unit arrives ready for use with the horizontal rotor, tube holders, and 2-year warranty included in the affordable price.

Tubes

6 75 & 100 mm 3-10 mL

Applications

Chemistry, Coag, PPP

Centrifugation

Horizontal 1,600 xg Max RCF 3,400 Max RPM

Included

Rotor Tube Holders


Warranty

2 years
Rotor: lifetime warranty
Extended warranty available

DRUCKER MODEL 642E

Part Numbers	642E
Centrifuge	00-018-009-001
75-100 mm tube holder	7713079 (1 qty)
SmartView Platform	00-079-009-001
1 yr Warranty Extension	03-0-0007-0001
2 yr Warranty Extension	03-0-0007-0002
3 yr Warranty Extension	03-0-0007-0006

Specifications	642E
Sound level	64 dBA
Dimensions	11.75 x 14 x 8.75 in 30 x 35.5 x 22 cm
Weight	11 lbs (5 kg)
Max Acceleration	10 seconds
Max Deceleration	20 seconds
Power	115 V – 50/60 Hz



OTHER DRUCKER MODELS



For mobile centrifugation For use with 12V outlet



HORIZON

3 preset cycles, available for 6, 12, 24 tubes For blood and urine



SampleSafe[™] Lockbox

24-hour sample protection Top load and front load models







RoHS Compliant



ISO 13485 certified



druckerdiagnostics.com sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only) +1-814-692-7661

