

PCDHGC3 Polyclonal Antibody

Catalog number: 10857-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:


PCDHGC3; PC 43, PC43, PCDH

GAMMA C3, PCDH2,

PCDHGC3, Protocadherin 2,

Protocadherin 43,

Protocadherin gamma C3



Immunohistochemical of paraffin-embedded human breast cancer using 10857-1-AP (PCDHGC3 antibody) at dilution of 1:100 (under 10x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PCDHGC3 MW:	150 kDa
Observed PCDHGC3 MW:	
Positive IHC detected in	Human breast cancer tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1243
GenBank accession number:	BC019299
Gene ID (NCBI):	5098
Full name:	Protocadherin gamma subfamily C, 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.