

**EDGE BIOLOGICALS, Inc.**

598 North Second Street  
Memphis, TN 38105  
Phone: 901-523-0034  
Fax: 901-527-3343

# Certificate of Analysis

**QUALITY ASSURANCE STATEMENT**

Quality control testing has been performed on each lot of culture medium according to Document M22-A3 Approved Standard, "Quality Assurance for Commercially Prepared Microbiological Culture Media" (CLSI) and internal EBI specifications.

<b>PRODUCT: Sterile Saline 0.85% 9.9 mL</b>		<b>LOT: 355</b>	
DATE PREPARED: 12/21/2018		EXP: 12/21/2019	
pH : 7.15	APPEARANCE: Clear	STERILIZATION: 121°C, 25 minutes	
Prepared from: DI Water	Lot: N/A	Exp: N/A	
Sodium Chloride (Fisher)	Lot: 8GFR118120033	Exp: N/A	

This medium can be adjusted or supplemented according to the performance criteria required.

**Quality Control:**

Sample of this product was enriched into tryptic soy broth and incubated for three days at 35 ± 2°C; result was No Growth. Aerobic Plate Count performed on a sample of this product and incubated for three days at 35 ± 2°C; result was 0 CFUs.

**Release Date: 12/27/2018**

Paul Horton  
Quality Control Manager

*Paul Horton*