

## EIF3L Polyclonal Antibody

Catalog number: 11482-1-AP

Size: 60 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

EIF3L; EIF3EIP, EIF3L, EIF3S11,

EIF3S6IP, HSPC021, HSPC025,

MSTP005



mouse brain tissue were subjected to SDS PAGE followed by western blot with 11482-1-AP( EIF3L Antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated EIF3L MW:	69 kDa
Observed EIF3L MW:	67 kDa
Positive WB detected in	Mouse brain tissue, rat brain tissue
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag1974
GenBank accession number:	BC007510
Gene ID (NCBI):	51386
Full name:	Eukaryotic translation initiation factor 3, subunit L

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.