

Product(s): Reveal® G4 Rapid HIV-1/2 Antibody Test

Subject: Maximum Expiry Date Extension

Dear Customers and Partners,

We are pleased to inform you that MedMira's Reveal® G4 Rapid HIV-1/2 Antibody Test expiry date will be increasing to 18 months from 12 months. There is no change of testing procedure steps nor interpretation of the results. An overview of the product specifications provided below.

Product Features:

- 18 month shelf life Simple 3-step procedure
- Temperature limit (2-30°C) Whole blood, serum or plasma specimens
- Instant results On-board procedural and reagent control line

Sensitivity - HIV-1

- Serum 99.8% Whole Blood (Fingerstick) 99.64%
- Plasma 99.8% Whole Blood (Venipuncture) 99.46%

Sensitivity - HIV-2

The overall sensitivity of the Reveal® G4 Rapid HIV-1/2 Antibody Test to detect HIV-2 Antibodies was calculated to be 100% (95% CL = 98.20% to 100.00%).

Specificity

- Serum 99.1% Whole Blood (Fingerstick) 99.71%
- Plasma 98.6% Whole Blood (Venipuncture) 99.90%

Packaging Formats

Product Format	Cat. No.	Specimen Type	Qty. per Box
Reveal® G4 POC	CHM06Z604	Fingerstick Whole Blood	20 tests
Reveal® G4 LAB+	CHM06Z603	Venipuncture Whole Blood Serum and Plasma	50 tests
Reveal® G4 LAB S/P	CHM06Z602	Serum and Plasma	50 tests

Should you have any questions regarding this Customer Advisory Notice or Reveal® G4 please contact our Customer Support Team at G4@medmira.com or +1 877 633 6472