

## SLC27A3 Polyclonal Antibody

Catalog number: 12943-1-AP


Size: 35 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

SLC27A3; ACSVL3, FATP 3,  
FATP3, Fatty acid transport  
protein 3, SLC27A3,  
UNQ367/PRO703, VLCS 3



human testis tissue were subjected to SDS PAGE followed by western blot with 12943-1-AP(SLC27A3 antibody) at dilution of 1:500



IP Result of anti-SLC27A3 (IP:12943-1-AP, 4µg; Detection:12943-1-AP 1:500) with mouse testis tissue lysate 4000µg.

### Applications

Tested applications:	ELISA, WB, IP
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated SLC27A3 MW:	75 kDa
Observed SLC27A3 MW:	79 kDa
Positive WB detected in	Human testis tissue, mouse testis tissue
Positive IP detected in	Mouse testis tissue
Recommended dilution:	WB: 1:1000-1:10000 IP: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4015
GenBank accession number:	BC029792
Gene ID (NCBI):	11000
Full name:	Solute carrier family 27 (fatty acid transporter), member 3

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.