

GEB Universal tips for TWC thin wallcrown And for improved flexibility.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (7.0*4.0)mm	Packaging Spec.
PT1000-B-B	1000ul	Blue				1000tips/bag, 5bags/case, 5000tips/case.
PT1000-B-B-LB	1000ul	Blue	Low Retention			
FT1000-B-N	1000ul	Natural			Filtered	
FT1000-B-N-LB	1000ul	Natural	Low Retention		Filtered	
Case dimension L*W*H: 690*235*330mm. Gross Weight: 5.8kg.						
PT1000-1-B	1000ul	Blue				96tips/rack, 10racks/pack, 5packs/case, 4800tips/case.
PT1000-1-B-LB	1000ul	Blue	Low Retention			
PT1000-1-BS	1000ul	Blue		Sterilized		
PT1000-1-BS-LB	1000ul	Blue	Low Retention	Sterilized		
FT1000-R-NS	1000ul	Natural		Sterilized	Filtered	
FT1000-R-NS-LB	1000ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 595*285*475mm. Gross Weight: 15kg.						

- ✧ 1000μL with graduation marks at 100ul ,200ul and 500ul
- ✧ 96 Tip/Hanging Tip Rack format / Bulk
- ✧ Polypropylene
- ✧ Presterile/Nonsterile
- ✧ Filtered (With aerosol-resistant filters) /NonFiltered
- ✧ Certified RNase-/DNase-free
- ✧ Nonpyrogenic
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Universal tips workstations.

Drawing:

