

FABP1 Polyclonal Antibody

Catalog number: 13626-1-AP


Size: 29 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

FABP1; FABPL; L-FABP



human heart tissue were subjected to SDS PAGE followed by western blot with 13626-1-AP(FABP1 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human heart tissue
Calculated FABP1 MW:	127aa,14kd
Observed FABP1 MW:	14kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4540
GenBank accession number:	BC032801
Gene ID (NCBI):	2168
Full name:	Fatty acid binding protein 1, liver

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.