

Langerin Polyclonal Antibody

Catalog number: 11841-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CD207; CD207, CD207 molecule,
langerin, CLEC4K, Langerin



human skin tissue were subjected to SDS PAGE followed by western blot with 11841-1-AP(CD207 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated Langerin MW:	328aa,37kd
Observed Langerin MW:	40kd
Positive WB detected in	Human skin tissue
Recommended dilution:	WB: 1:250 - 1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2479
GenBank accession number:	BC022278
Gene ID (NCBI):	50489
Full name:	CD207 molecule, langerin

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.