

GLS Aqua Plus Series Linear Shaking Baths

Grant

**Precision engineered and
newly designed for 2010**

**Precision
temperature controls
with convenient
temperature and
shaking presets**



- Robust magnetically coupled shaking mechanism – safe, reliable, durable, and space saving
- Drain tap – allows convenient tank emptying
- Two programmable temperature and shaking presets – allows commonly used operation parameters to be quickly selected
- Countdown timer with audible buzzer – enables accurate reaction timing
- Grant non-drip polycarbonate lid included as standard
- Choice of two models – 12 and 18 liter
- Three year warranty as standard

Model number	Tank size liters
GLSAQP12US	12
GLSAQP18US	18

**For complete features and specifications
visit our website at www.grantsci.com**

Grant Instruments
601 Route 206, Suite 26-730
Hillsborough, NJ 08844
800-992-3028 • Fax: 908-431-3029

SBB Aqua Plus Series Boiling Baths

Grant

**Precision engineered and
newly designed for 2010**

**New robust design
for continuous
operation at 100°C**



- Stainless steel tank in a robust outer case – tough and durable in demanding environments
- Non-drip polycarbonate lid included as standard – improves performance, limits evaporation, and conserves energy
- Constant level device maintains liquid level – ensures continuous optimal operation
- Two over temperature cutouts – protects users and workplace if bath is accidentally run without liquid
- Choice of 4 sizes to meet individual application requirements
- Flat lids with ring sets available as optional accessory
- Three year warranty as standard

Model number	Tank size liters
SBBAQP5US	5
SBBAQP12US	12
SBBAQP18US	18
SBBAQP26US	26

**For complete features and specifications
visit our website at www.grantsci.com**

Grant Instruments
601 Route 206, Suite 26-730
Hillsborough, NJ 08844
800-992-3028 • Fax: 908-431-3029