

## MYO1D Polyclonal Antibody

Catalog number: 14028-1-AP


Size: 45 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MYO1D; KIAA0727, MYO1D,  
myosin ID, myr4



mouse lung tissue were subjected to SDS PAGE followed by western blot with 14028-1-AP(MYO1D antibody) at dilution of 1:200

### Applications

Tested applications:	ELISA, WB
Cited applications:	IHC
Species specificity:	Human, Mouse; other species not tested.
Cited species:	Human
Calculated MYO1D MW:	116 kDa
Observed MYO1D MW:	100kd
Positive WB detected in	Mouse lung tissue, mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5176
GenBank accession number:	BC030602
Gene ID (NCBI):	4642
Full name:	Myosin ID

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.