

PIP5K3 Polyclonal Antibody

Catalog number: 13361-1-AP

Size: 62 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

PIP5K3; CFD; FAB1; KIAA0981;

MGC40423; PIKFYVE; PIP5K



HeLa cells were subjected to SDS PAGE followed by western blot with 13361-1-AP(PIP5K3 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HeLa cells, mouse brain tissue, mouse lung tissue
Calculated PIP5K3 MW:	451aa,50kd
Observed PIP5K3 MW:	237kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4177
GenBank accession number:	BC032389
Gene ID (NCBI):	200576
Full name:	Phosphatidylinositol-3-phosphate/phosphatidylinositol 5-kinase, type III

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.