

Diamond® White Glass microscope slides




COLOR CODED SLIDES

Diamond® White Glass color coded microscope slides are offered in bright, attractive colors. They are available with ground edges and 90° corners, or with beveled edges and clipped corners. All are printed with "Globe" under the color area for orientation purposes. They are pre-cleaned and ready to use.

- Color coded on one end, one side
- Produced from the finest quality optical white glass
- Packed in high-quality plastic boxes
- Packaging: 72/plastic box, 20 boxes/case (10 gross)

Ground Edge with 90° Corners

Thomas #	Item #	Description	Unit
1208T55	1380-50B	Frosted, blue	1440
1208T57	1380-50G	Frosted, green	1440
1208T58	1380-50L	Frosted, lilac	1440
1208T60	1380-50P	Frosted, pink	1440
1208T62	1380-50W	Frosted, white	1440
1208T63	1380-50Y	Frosted, yellow	1440



CHARGED SLIDES

Globe's charged microscope slides provide the ideal surface for cell and tissue adhesion. These popular slides are ideal for use in histology, cytology and pathology departments, where the convenience of an adhesive slide saves the user time and money. They are pre-cleaned and ready to use.

- Adhesive surface
- Color coded on one end, one side
- Packaging: 72/plastic box, 20 boxes/case (10 gross)

Ground Edge with 90° Corners

Thomas #	Item #	Color	Unit
1208T66	1358A	Frosted, aqua	1440
1208T67	1358B	Frosted, blue	1440
1207N24	1358D	Frosted gold	1440
1208T68	1358G	Frosted, green	1440
1208T69	1358L	Frosted, lilac	1440
1208T70	1358N	Frosted, orange	1440
1208T71	1358P	Frosted, pink	1440
1208T72	1358T	Frosted, tan	1440
1207N23	1358W	Frosted, white	1440
1208T73	1358Y	Frosted, yellow	1440

Beveled Edge with Clipped Corners

Thomas #	Item #	Color	Unit
1186N61	1354W	Frosted, white	1440

BUY 3 CASES of Diamond White Glass Microscope Slides between April 1 and June 30, 2025, and GET 1 CASE FREE!


Limit 3 free cases. Mix & match. Free item of equal or lesser value. Offer good while supplies last. Valid for shipment within the contiguous US only.

For FREE Product(s), Please Complete Redemption Information on Reverse Side!

**Buy 3
Get 1
FREE**



COLOR CODED SLIDES



CHARGED SLIDES

ACCOUNT INFORMATION

First Name Last Name Date

Organization Account No.

Shipping Address Rm./Bldg./Dept.

City State Zip Code

Phone Email



**MY FREE Product
Selection(s)**

Allow 3-4 weeks for delivery of
free products.

Valid 4/01/25 - 6/30/25

Complete the account information and
submit, along with Thomas Scientific proof of
purchase, to:

Promotions@thomassci.com

Reference **DWHITE3/1** when submitting redemption.

*Disclaimer: Limit 3 free cases. Mix & match. Free item of equal or
lesser value. Offer good while supplies last. Valid for shipment
within the contiguous US only. Redemptions must be received
no later than July 15, 2025. Pricing and offers subject to change
without notice.*