

## PCDHB12 Polyclonal Antibody

Catalog number: 14020-1-AP

Size: 48 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PCDHB12; PCDH beta 12, PCDH

BETA12, PCDHB12,

protocadherin beta 12



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14020-1-AP(PCDHB12 antibody) at dilution of 1:800

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated PCDHB12 MW:	87 kDa
Observed PCDHB12 MW:	87kd
Positive WB detected in	Mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5135
GenBank accession number:	BC045637
Gene ID (NCBI):	56124
Full name:	Protocadherin beta 12

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.