

| | |
|---------------|--|
| Product Name: | ABS ET 8 |
| Description: | <p>ABS ET Series Includes 6 Canisters ½ CC Straw (1 Level Bulk): 1122 Static Evaporation Rate (Liters/day): 0.15 LN2 Capacity: 8.4 Normal Working Time: 35 Days Neck Diameter (mm): 55.4 Exterior Height (mm): 469.9 Exterior Diameter (mm): 259.1 Canister Height (mm): 127 Canister Diameter (mm): 41.9 Designed for small capacity storage Durable aluminum tank with powder coated exterior High strength neck tube reduces liquid nitrogen loss Lid lock provision Low liquid nitrogen consumption Liquid nitrogen measuring rod for level validation One-year warranty on all parts and labor Three-year vacuum warranty</p> |