

## MBD4 Polyclonal Antibody

Catalog number: 11270-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MBD4; MED1



human skin tissue were subjected to SDS PAGE followed by western blot with 11270-1-AP(MBD4 Antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse; other species not tested.
Positive WB detected in:	Human skin tissue, A549 cells
Calculated MBD4 MW:	66kd
Observed MBD4 MW:	66kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag1822
GenBank accession number:	BC011752
Gene ID (NCBI):	8930
Full name:	Methyl-CpG binding domain protein 4

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.