

LCAT Polyclonal Antibody

Catalog number: 12243-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

LCAT; LCAT



mouse testis tissue were subjected to SDS PAGE followed by western blot with 12243-1-AP(LCAT antibody) at dilution of 1:500

Applications

| | |
|-------------------------|---|
| Tested applications: | ELISA, WB |
| Species specificity: | Human, mouse, rat; other species not tested. |
| Calculated LCAT MW: | 440aa,50 kDa |
| Observed LCAT MW: | 55kd |
| Positive WB detected in | Mouse testis tissue, mouse bladder tissue, PC-3 cells |
| Recommended dilution: | WB: 1:500-1:5000 |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

| | |
|---------------------------|--------------------------------------|
| Immunogen: | Ag2887 |
| GenBank accession number: | BC014781 |
| Gene ID (NCBI): | 3931 |
| Full name: | Lecithin-cholesterol acyltransferase |

Product information

| | |
|----------------------|--|
| Purification method: | Antigen affinity purification |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |