




Specifications

Electrical requirements	110VAC 60Hz single phase, 500 watts
Chamber dimensions	10 x 9 x 10" (WxDxH)
Unit dimensions	21 x 17 x 17" (WxDxH)
Shipping dimensions	26 x 23 x 22" (WxDxH)
Temperature range	Ambient temperature to 400°F
Temperature uniformity	+/- 2°F
Temperature controller precision	0.1°F
Heating elements	In the back of the oven
Timer range	1 to 9999 minutes
Racks included	2 (additional racks available at extra cost)
Maximum rack capacity	3 (height between each rack is 2 inches)
Exhaust outlet diameter	1 1/4"
Unit weight	55 Lb
Shipping weight	100 Lb
Access port with cover	Optional 1 1/2 or 2" diameter access port at \$290
Certification	CE, UL/CSA (optional)

Standard package

Part description	Quantity	Part image
AI 0.5 cu ft forced air oven	1 pc	
Stainless steel rack and brackets	2 sets	
110 volt 3-prong power cord	1 pc	
User's manual	1 pc	