

FAM3D Polyclonal Antibody

Catalog number: 12336-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

FAM3D; FAM3D, OIT1, Protein

FAM3D, EF7



mouse colon tissue were subjected to SDS PAGE followed by western blot with 12336-1-AP(FAM3D antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated FAM3D MW:	224aa,24 kDa
Observed FAM3D MW:	30kd
Positive WB detected in	Mouse colon tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2965
GenBank accession number:	BC015359
Gene ID (NCBI):	131177
Full name:	Family with sequence similarity 3, member D

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.