



P.O. Box 240746  
Charlotte, NC 28224  
www.reagents.com

info@reagents.com  
Phone: 704.554.7474  
Toll Free: 800.732.8484  
Fax: 888-843-4384

## CERTIFICATE OF ANALYSIS

Product: **Sodium Salicylate, Reagent Grade**

Date of Manufacture: 1/1/2018

Reference ID: C1324700

Expiration Date: January 2028

Product Code: C1324700

Recommended Storage: 15 - 30 C (59°F - 86°F)

Lot Number: **7905557**

CAS Number: 54-21-7

Test	Specification	Result
Assay	>= 99.0 %	99.5%
Heavy Metals (as Pb)	<= 0.002 %	0.001%
Water (H <sub>2</sub> O)	<= 0.5 %	0.1%

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Quality Control

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials \*\* Contents of Certificates and Labels."

2019052408:53:15sfarrell-0-0