



# Certificate of Analysis

Created May 27, 2022

**Product:** Human Serum Albumin

Powder

**Catalog #: 800-125P**

**Lot #: H122028**

**Origin: United States (100%)**

**Manufacture Date: May2022**

**Expiration Date: Sep2028**

TEST	METHODOLOGY	SPECIFICATION	RESULTS
<b>Microbiological Testing</b>			
Microbial Growth		<2.99 cfu/mL	0 cfu/mL
<b>Endotoxin Testing</b>			
Endotoxin		≤2 EU/mg	0.96 EU/mg
<b>Viral Testing</b>			
HBsAg		Non-Reactive	Non-Reactive
Anti-HCV		Non-Reactive	Non-Reactive
Anti-HIV-1/HIV-2		Non-Reactive	Non-Reactive
Syphilis		Negative	Negative
HIV-NAT		Not Detected	Not Detected
HCV-NAT		Not Detected	Not Detected
<b>Physical Testing</b>			
pH @ RT			7.1
<b>Physical Testing</b>			
Purity			97.0%

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use. Results shown were obtained by carefully performed methods believed to be reliable. However, since some results may vary for specific tests depending upon methodology and other variables, it is suggested that tests for which results are particularly important be repeated by the user of this product.*

*All human plasma is collected from stringently screened donors at FDA-licensed collection centers in the United States.*

VANESSA YOUNG 27MAY2022

Created By/ date

Qu 27MAY 2022

Verified By/ date