

TBRG1 Polyclonal Antibody

Catalog number: 13346-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TBRG1; FLJ14621; FLJ25020;

FLJ90113; MGC129890; NIAM;

TB-5



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 13346-1-AP (TBRG1 Antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	BxPC-3 cells
Calculated TBRG1 MW:	45kd
Observed TBRG1 MW:	45kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4179
GenBank accession number:	BC032312
Gene ID (NCBI):	84897
Full name:	Transforming growth factor beta regulator 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.