

Features:

- ✓ Large LCD display provides all visual information;
- ✓ Unique temperature display design for easily read;
- ✓ Whisper quiet operation, suit for lab environment;
- ✓ Brushless motor drive, maintenance free;
- ✓ Affordably priced, fits into your budget;
- ✓ Bio-containment rotor, safe and autoclavable;
- ✓ Special double cooling air design, lower rotor heating (D3024);
- ✓ Maintains 4°C at maximum speed to protect precious samples (D3024R);
- ✓ Powerful compressor and cooling system, fast cooling (D3024R).
- ✓ Additional rotors include a 4 x 8-tube PCR strip rotor



D3024 Micro Centrifuge



D3024R Refrigerated Micro Centrifuge

Specifications

Max. Speed	15000rpm (300-15000rpm), Increment: 10rpm
Max. RCF	21380×g, Increment: 10×g
Rotor	2ml×24, 0.5ml×36, PCR8 serial tube
Run Time	1min~99min Timer or Continuous Operation
Temperature Range	-16°C~40°C (D3024R)
Refrigerant	HFC134a
Driving Motor	Brushless Motor
Safety Devices	Door Interlock, Over-Speed Detector, Over-Temperature Detector, Automatic Internal Diagnosis
Power	Single-Phase, ~110V, 60Hz, 5A
Dimensions (inches)	(W) 9.4× (L) 13× (H) 9.8 (D3024) (W) 11.4× (L) 20.4× (H) 11.7 (D3024R)
Weight	66lbs (D3024R) 27lbs (D3024)
Additional Features	Speed/Acceleration switch Short-Time Run Function Sound-alert Function

2 YEAR WARRANTY

Ordering information

Cat No.	Description
912000069999	SCIOGEX D3024 High Speed Micro-Centrifuge, c/w AS24-2 24-Place Rotor, 110V/60Hz
922000069999	SCIOGEX D3024R High Speed Refrigerated Micro-Centrifuge, c/w AS24-2 24-Place Rotor, 110V/60Hz
19400002	24 Place Rotor kit with cover, fixing clips and O' ring, 15000rpm, 2ml x 24, used with D3024 & D3024R
19400003	36 Place Rotor kit with cover, fixing clips and O' ring, 15000rpm, 0.5ml x 36, used with D3024 & D3024R
19400004	4-PCR8 Rotor kit with cover, fixing clips and O' ring, 15000rpm, 4-PCR8, used with D3024 & D3024R
19500001	0.2ml rotor adapter, used with 19400002 24 place rotor, 24pcs/pk
19500002	0.5ml rotor adapter, used with 19400002 24 place rotor, 24pcs/pk
19500003	0.2ml rotor adapter, used with 19400003 36 place rotor, 36pcs/pk

Accessories



Rotor Cat. No.		Max. Speed (rpm)	Max. Rcf (×g)	Capacity	Bio-Safe	Centrifuge	
						D3024	D3024R
19400002	Base info.	15000	21380	2ml×24	Yes	●	●
	Adapters	19500001 - 0.2ml 19500002 - 0.5ml					
19400003	Base info.	15000	21380	0.5ml×36	Yes	●	●
	Adapters	19500003 - 0.2ml					
19400004	Base info.	15000	21380	4-PCR8	Yes	●	●