

Corning® MSCulture Max™ Media

CORNING

Over the past decade, stem cell research has made significant strides and has increasingly been applied as a therapeutic method for various disorders. Stem cells, by definition, are specialized cells naturally found in most body tissues, capable of self-renewal and differentiation into diverse cell types.

Various types of stem cells are under investigation for therapeutic applications, and their progress towards clinical use varies. Among these, mesenchymal stem/stromal cells (MSCs) have been leading the way, outpacing other multipotent and pluripotent cell types. MSCs boast a diverse biological toolkit that can be harnessed for therapeutic purposes.

Corning's new MSCulture Max media increase final hMSC yields, streamlines culture, reduces in-process risk, and eliminates media exchanges, enabling higher cell expansion with less media and fewer media changes. The MSCulture Max media is available in 3 versions.

- ▶ Standard research-grade medium for MSC expansion
- ▶ Xeno-free research-grade medium optimized for rapid MSC expansion
- ▶ GMP-grade medium with pathogen-reduced components for clinical application


NOTE: To achieve optimal MSC expansion performance, you will need to order both the basal medium and the corresponding supplement individually, which are outlined in the Ordering Information section. After ordering the basal medium and supplement separately, both components must be combined to make MSCulture Max media complete per the Instructions for Use (CLS-AN-864DOC).

Features

- ▶ High-performance expansion media for batch culture systems
- ▶ No media exchanges required
- ▶ No surface coating required

Benefits

- ▶ **Unparalleled productivity:** Achieves over 10-fold hMSC expansion within 4-5 days without media changes, significantly reducing time and labor costs.
- ▶ **Transition-friendly formulation:** Research, xeno-free and GMP versions available for rapid translation pathway from research into clinical manufacturing.
- ▶ **Scalable and efficient bioprocessing:** Enables high-volume cell production with minimal media consumption.



Ordering Information

Corning® MSCulture Max™ Media

NOTE: To achieve optimal MSC expansion performance, you will need to order both the basal medium and the corresponding supplement individually.*

Cat. No.	Description	Shelf Life of Mixed Basal Medium and Corresponding Supplement*	Shelf Life	Shipping Method**
Corning MSCulture Max Media, Standard Research-grade				
42-010-CV	MSCulture Max basal medium	2 weeks	15 months	Refrigerated
42-000-CR	MSCulture Max supplement		2 years	Frozen
Corning MSCulture Max-XF Media, Xeno-free Research-grade				
42-010-CV	MSCulture Max basal medium	2 weeks	15 months	Refrigerated
42-100-CR	MSCulture Max-XF supplement		2 years	Frozen
Corning MSCulture Max-GMP Media, GMP-grade for Bioproduction				
42-010-CV	MSCulture Max basal medium	1 week	15 months	Refrigerated
42-200-CR	MSCulture Max-GMP supplement		3 years	Frozen

*As stated above, both components must be combined to make MSCulture Max media complete per the Instructions for Use (CLS-AN-864DOC) to achieve optimal MSC expansion performance.

**You will receive two different shipments arriving separately, which may occur on different dates.

How to Purchase: For specific availability in your region and purchasing options, terms and conditions of sale, customer/product support, and certificates, visit www.corning.com/how-to-buy.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.



Corning Incorporated
Life Sciences
www.corning.com/lifesciences

NORTH AMERICA
t 800.492.1110
t 978.442.2200

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India
t 91 124 4604000
Japan
t 81 3-3586 1996
Korea
t 82 2-796-9500
Singapore
t 65 6572-9740
Taiwan
t 886 2-2716-0338

EUROPE
CSEurope@corning.com
France
t 0800 916 882
Germany
t 0800 101 1153
The Netherlands
t 020 655 79 28
United Kingdom
t 0800 376 8660

All Other European Countries
t +31 (0) 206 59 60 51
LATIN AMERICA
grupolA@corning.com
Brazil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400