



Note: labeling is forbidden on the Freezing Bags body;
a small label can be placed on the surface of the inner PE bag
and a product label on the outer PE bag

Item No.	Qty	Description	Inventory Code	Material
D1	1	50mL SU Freezing Bags (0.3)	FB28EVA120	EVA
A1	2	0.203"ID/0.280"OD(5.15mm/7.1mm)EVA	YP001	EVA
B1	1	2.95*4.15mm EVA Tubing,Length 18cm	JNTH4030	EVA
C1	1	Freeze tube tubing assembly	JNSD002	PVC/Silicone/PE
F1	1	Light-curing glue	1-CN007	PAA
F2	2	PE Bag	WJM1001	PE
F3	1	2# Carton,,5pcs/pack, 20 packs/box	BZX01A	Corrugated Paper

Customer Approval		LePure Biotech Co., Ltd.			
Company:		Drawing No. CR0050-3E			
Signature:		Title: SU Freezing Bags			
Date:		Designed By:	Date:	Checked By:	Date: Edition
		Carl	2023.05.11	Zhou	2023.05.11 Rev01
Design Approval		We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property			
the manufacturing of the system(s) based on this design. By signing, you also assume full					