


Microplate Seals

Film and Foil Microplate Heat Seals

A comprehensive range of film and foil microplate heat seals, ideal for use with the Vitl VTS and Vitl Micro^{TS}



Description

Vitl stock a comprehensive range of film and foil microplate heat seals, ideal for (but not limited to) use with the Vitl VTS Microplate Heat Sealer and Vitl Micro^{TS} Microplate Heat Sealer.

All Vitl microplate seals are supplied in packs of 100 sheets, each measuring 125mm x 80mm.

Reels of seals are available on request.

Technical Specifications

	Peelable?	Pierceable?	Temperature range	Seals Polypropylene?	Seals Polyethylene?	Seals Polystyrene?	Seals COC?	Color code	Catalogue Number	Product Notes
20µm Heat Sealing Foil	X	✓	-20 - +120°C	✓	✓	✓	X	Green	V901002	
38µm Heat Sealing Foil	✓	✓	-20 - +110°C	✓	✓	✓	X	Blue	V901003	Only peels from PS
70µm Heat Sealing Foil	✓	X	-80 - +80°C	✓	✓	X	X	Red	V901001	DMSO resistant
75µm Clear Sealing Film	X	✓	-80 - +80°C	✓	✓	✓	✓	75µm	V901005	Suitable for PCR at 110°C
85µm Heat Sealing Foil	✓	X	-200 - +110°C	✓	X	X	X	Black	V901004	
105µm Clear Sealing Film	✓	X	-80 - +180°C	✓	✓	✓	✓	105µm	V901006	Suitable for PCR at 110°C

Related Products



VTS
Variable
Temperature
Sealer
V902001



Micro^{TS}
Microplate
Heat Sealer
V903001

Vitl 
LIFE SCIENCE SOLUTIONS
A member of the ITL group