

## COX6A2 Polyclonal Antibody

Catalog number: 11421-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

COX6A2; COX VIa M, COX6A,

COX6A2, COX6AH, COXVIAH

### Background

COX6A2(Cytochrome c oxidase subunit 6A2, mitochondrial) is also named as COX6A, COX6AH and Belongs to the cytochrome c oxidase subunit 6A family. This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

### Applications

**Tested applications:** ELISA  
**Cited applications:** WB  
**Species specificity:** Human,Mouse,Rat; other species not tested.  
**Cited species:** Human  
**Calculated COX6A2 MW:** 97aa,11 kDa  
**Observed COX6A2 MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

**Immunogen:** Ag1985  
**GenBank accession number:** BC029818  
**Gene ID (NCBI):** 1339  
**Full name:** Cytochrome c oxidase subunit VIa polypeptide 2

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.