

# METTLER TOLEDO



Karl Fischer Titrators	V10S	V20S	V30S	C10SD/C10SX	C20SD/C20SX	C30SD/C30SX
<b>Catalog Number</b>	30267116	30252666	30252667	30267115 (SD) / 30267114 (SX)	30252661 (SD) / 30252660 (SX)	30252663 (SD) / 30252662 (SX)
<b>Hardware</b>						
<b>Burette Sizes (mL)</b>	1, 5, 10, 20	1, 5, 10, 20	1, 5, 10, 20	N/A	N/A	N/A
<b>Burette Drive(s)</b>	1	1	1	N/A	N/A	N/A
<b>Standard Delivery Burette</b>	DV1005 5mL	DV1005 5mL	DV1005 5mL	N/A	N/A	N/A
<b>Standard Delivery Electrode(s)</b>	DM143-SC	DM143-SC	DM143-SC	Generator electrode & DM143-SC	Generator electrode & DM143-SC	Generator electrode & DM143-SC
<b>Automation Option(s)</b>	✗	✗	InMotion KF Flex/Pro Oven Autosamplers	✗	✗	InMotion KF Flex/Pro Oven Autosamplers
<b>Solvent Manager</b>	✓	✓	✓	Optional	✓	✓
<b>User Interface</b>						
<b>Titration Types</b>	Karl Fischer Volumetric	Karl Fischer Volumetric	Karl Fischer Volumetric	Karl Fischer Coulometric	Karl Fischer Coulometric	Karl Fischer Coulometric & Bromine Index Determination
<b>Max Number of User Methods</b>	5	5	150	5	5	150
<b>Shortcuts Per User</b>	4	12	12	4	12	12
<b>Connectivity</b>						
<b>Data Export</b>	USB, Network, PDF	USB, Network, PDF	USB, Network, PDF	USB, Network, PDF	USB, Network, PDF	USB, Network, PDF
<b>Printers</b>	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn
<b>PC Software/Connectivity</b>	N/A	LabX Express or LabX Server	LabX Express or LabX Server	N/A	LabX Express or LabX Server	LabX Express or LabX Server
<b>Compatible Smart Accessories</b>	SmartChemical Reader	SmartChemical Reader	SmartChemical Reader	N/A	N/A	N/A
<b>Services</b>						
<b>EQPac + SCC Calibration</b>	✓	✓	✓	✓	✓	✓
<b>IPac (Calibration Optional)</b>	✓	✓	✓	✓	✓	✓
<b>VPac</b>	✗	✗	✗	✗	✗	✗
<b>Extended Care (Calibration Optional)</b>	✓	✓	✓	✓	✓	✓
<b>On-Site Service</b>	✓	✓	✓	✓	✓	✓

# METTLER TOLEDO

## Compact Titrators



**G10S**



**G20S**

## Titration Excellence



**T5**



**T7**



**T9**

Catalog Number	30267117	30252669	30252672	30252675	30252676
<b>Hardware</b>					
Burette Sizes (mL)	1, 5, 10, 20	1, 5, 10, 20	1, 5, 10, 20	1, 5, 10, 20	1, 5, 10, 20
Burette Drive(s) to Dose & Titrate	1	1 internal + 1 external - dosing only	1 internal + 1 external	1 internal + 3 external	1 internal + 7 external
Standard Delivery Burette	DV1020 20mL	DV1020 20mL	DV1010 10 mL	Not Included	Not Included
Standard Delivery Electrode(s)	DGi115-SC	DGi115-SC	Not included	Not Included	Not Included
Automation Option(s)	Rondolino	Rondolino	Rondolino, InMotion Autosamplers	Rondolino, InMotion Autosamplers, InMotion KF Flex/Pro autosamplers	Rondolino, InMotion Autosamplers, InMotion KF Flex/Pro autosamplers
Solvent Manager	Optional	Optional	Optional	Optional	Optional
<b>User Interface</b>					
Titration Types	Acid/Base, Redox, Complexometric, Precipitation, Photometric	Acid/Base, Redox, Complexometric, Precipitation, Photometric	Acid/Base, Amperometric, Bromine Index Determination, Complexometric, Conductivity, Karl Fischer Coulometric, Karl Fischer Volumetric, Precipitation, Photometric, Redox, Standard Addition, Thermometric, Voltametric	Acid/Base, Amperometric, Bromine Index Determination, Complexometric, Conductivity, Karl Fischer Coulometric, Karl Fischer Volumetric, Precipitation, Photometric, Redox, Standard Addition, Thermometric, Voltametric	Acid/Base, Amperometric, Bromine Index Determination, Complexometric, Conductivity, Karl Fischer Coulometric, Karl Fischer Volumetric, Precipitation, Photometric, Redox, Standard Addition, Thermometric, Voltametric
Max Number of User Methods	5	150	150	150	150
Shortcuts Per User	4	12	24	24	24
<b>Connectivity</b>					
Data Export	USB, Network, PDF	USB, Network, PDF	USB, Network, PDF	USB, Network, PDF	USB, Network, PDFa
Printers	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn	USB-P25, OKi B412dn
PC Software/Connectivity	N/A	LabX Express	LabX Express or LabX Server	LabX Express or LabX Server	LabX Express or LabX Server
Compatible Smart Accessories	SmartChemical Reader	SmartChemical Reader	SmartChemical Reader, SmartSample	SmartChemical Reader, SmartSample	SmartChemical Reader, SmartSample
<b>Services</b>					
EQPac + SCC Calibration	✗	✗	✓	✓	✓
IPac (Calibration Optional)	✓	✓	✓	✓	✓
VPac	✗	✗	✗	✗	✗
Extended Care (Calibration Optional)	✓	✓	✓	✓	✓
On-Site Service	✓	✓	✓	✓	✓