

METTLER TOLEDO



EasyPlus Titrators	Easy pH	Easy Cl	Easy Ox	Easy Pro	Easy KFV	Easy EPM	Easy Na Sodium Analyzer
Catalog Number	30060041	30060043	30060042	30060044	30060045	30470855	30060051
Hardware							
Burette Sizes (mL)	10, 20, 50	10, 20, 50	10, 20, 50	10, 20, 50	10, 20, 50	10, 20, 50	10, 20, 50
Burette Drive(s)	1	1	1	1	1	1	1
Standard Delivery Burette	EasyPlus Burette 20mL	EasyPlus Burette 20mL	EasyPlus Burette 20mL	EasyPlus Burette 20mL	EasyPlus Burette 10mL	EasyPlus Burette 20mL	EasyPlus Burette 10mL
Standard Delivery Electrode(s)	EG11-BNC	EM45-BNC	EM40-BNC	EG11-BNC	EM43-BNC	Not Included	DX222-Na & DX205-SC
Automation Option(s)	×	×	×	×	×	×	×
Solvent Manager	×	×	×	×	×	×	×
User Interface							
Titration Types	Acid/Base	Precipitation	Redox	Acid/Base, Precipitation, Redox	Karl Fischer Volumetric	Manual Titrations with color indicator or mV electrode	Sodium Determination by Standard Addition
Max Number of User Methods	1	1	1	3	1	1	1
Shortcuts Per User	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Connectivity							
Data Export	PDF and CSV via USB Flash Drive	PDF and CSV via USB Flash Drive	PDF and CSV via USB Flash Drive	PDF and CSV via USB Flash Drive	PDF and CSV via USB Flash Drive	PDF and CSV via USB Flash Drive	PDF and CSV via USB Flash Drive
Printers	USB-P25	USB-P25	USB-P25	USB-P25	USB-P25	USB-P25	N/A
PC Software/Connectivity	EasyDirect	EasyDirect	EasyDirect	EasyDirect	EasyDirect	N/A	EasyDirect
Compatible Smart Accessories	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Services							
EQPac + SCC Calibration	×	×	×	×	×	×	×
IPac (Calibration Optional)	×	×	×	×	×	×	×
Extended Care (Calibration Optional)	×	×	×	×	×	×	×
On-Site Service	×	×	×	×	×	×	×
VPac	✓	✓	×	✓	×	×	×