

Product Name:	PH-ABT-HC-UCBI-0204-LH
Description:	<p>Premier Pharmacy/Vaccine Undercounter Refrigerator Built-In</p> <p>2.5 Cu. Ft. Capacity</p> <p>One solid swing door, left hinged</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)</p> <p>Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable</p> <p>Audible and visual high and low temperature alarms</p> <p>2 probes; 1 in air and 1 in glycol bottle</p> <p>Remote alarm contacts</p> <p>Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall</p> <p>Back wall evaporator cover</p> <p>Vaccine storage power cord warning label</p> <p>Forced draft circulation</p> <p>Self-closing door</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>LED light</p> <p>Probe access port (3/8")</p> <p>Powder coat exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Pharmacy refrigerator/freezer toolkit and temperature logs</p> <p>Two year parts and labor warranty</p> <p>Pyxis®, Omnicell® and AcuDose RX® compatible</p> <p>Exterior dimensions: 17 3/4"W x 22 3/4"D x 30 3/4"H</p> <p>Energy Star Certified for Pharmacy Refrigerators and Freezers</p> <p>Cycle defrost</p> <p>115V, 60 Hz, .90 Amps, 1/8 HP</p> <p>Energy Consumption 0.58 kWh/day</p> <p>ADA compliant</p> <p>UL/C-UL listed</p> <p>Shipping weight: 85 lbs.</p>