

ENTPD3 Polyclonal Antibody

Catalog number: 13021-1-AP

Size: 49 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ENTPD3; CD39L3; FLJ93839;

HB6; NTPDase-3



human brain tissue were subjected to SDS PAGE followed by western blot with 13021-1-AP(ENTPD3 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human brain tissue, Jurkat cells, Raji cells
Calculated ENTPD3 MW:	529aa,59kd
Observed ENTPD3 MW:	59kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4114
GenBank accession number:	BC029869
Gene ID (NCBI):	956
Full name:	Ectonucleoside triphosphate diphosphohydrolase 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.