

Product Name:	ABS VS 6
Description:	<p> ABS VS Class Vapor Shipper Number of Canisters: 5 Total Vial Capacity (5/Cane): 840 Liquid Nitrogen Capacity: 6.2 Liters ½ CC Straws (1 Level Bulk): 3750 ½ CC straws (10/Cane): 2500 Static Evaporation Rate (Liters/day): 0.35 Static Holding Time (days): 16 Neck Diameter (mm): 96.7 Exterior Height (mm): 559 Exterior Diameter (mm): 368.3 Canister Height (mm): 279.4 Canister Diameter (mm): 81.3 Optional protective shipping container ensure upright shipping High strength neck tube reduces liquid nitrogen loss Lid with locking provision Hydrophobic liquid nitrogen absorbent system 2-hour charging capability One-year parts and labor warranty One-year vacuum warranty </p>